
GENETIC TEST:
Norrie disease (NDP gene)

FULL NAME:	Norrie disease (NDP gene)
TEST TYPE:	Clinical
TEST SPECIALTY:	Molecular Genetics
TEST PURPOSE:	Carrier diagnosis, Mutation confirmation, Post-natal Diagnosis, Pre-implantation genetic diagnosis, Prenatal diagnosis
SPECIMEN:	Peripheral (whole) blood on EDTA, Amniotic fluid, Chorionic villi, Cell culture
METHOD CATEGORY:	Sequence analysis: entire coding region Mutation screening and sequence analysis of selected exons
METHOD TECHNIQUE:	Bi-directional Sanger Sequence analysis
RIZIV CODE:	565456-565460

EQA:	• DNA Sequencing - Sanger, • DNA Sequencing - Sanger, • DNA Sequencing - Sanger, • DNA Sequencing - Sanger, • DNA Sequencing - Sanger, • DNA Sequencing - Sanger, • DNA Sequencing - Sanger , • DNA Sequencing - Sanger
TURNAROUND TIME (MAXIMUM):	3 months (10 working days for prenatal diagnosis)
CREATED:	27 Aug 2019 - 13:25
CHANGED:	09 Mar 2023 - 16:21
URL:	https://laboguide.uzbrussel.be/laboguide#Analyses:norrie&&&&1054&COLLECTION&0

Source URL: http://gentest.healthdata.be/genetic_test/765

RELATED CONTENT

Related Diseases

- [Norrie disease](#)

Related Laboratories

- [Centrum Medische Genetica - UZ Brussel VUB](#)

Related Analytes

- [norrin cystine knot growth factor NDP](#)

Source URL: http://gentest.healthdata.be/genetic_test/765