

## GENETIC TEST: Hemophilia A

|                         |                                                                                                       |
|-------------------------|-------------------------------------------------------------------------------------------------------|
| <b>FULL NAME:</b>       | Hemophilia A                                                                                          |
| <b>DESCRIPTION:</b>     | Screening of all coding exons of the F8 gene.                                                         |
| <b>TEST TYPE:</b>       | Clinical                                                                                              |
| <b>TEST SPECIALTY:</b>  | Molecular Genetics                                                                                    |
| <b>TEST PURPOSE:</b>    | Carrier diagnosis,<br>Post-natal Diagnosis,<br>Prenatal diagnosis                                     |
| <b>SPECIMEN:</b>        | Peripheral (whole) blood on EDTA,<br>Amniotic fluid,<br>Chorionic villi                               |
| <b>METHOD CATEGORY:</b> | Targeted variant analysis<br>Sequence analysis: entire coding region<br>Deletion/duplication analysis |

|                                   |                                                                                                                                                   |
|-----------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>METHOD TECHNIQUE:</b>          | PCR based technique<br>Next Generation Sequencing (NGS)<br>Bi-directional Sanger Sequence analysis<br>MLPA based techniques                       |
| <b>RIZIV CODE:</b>                | 565471-565482                                                                                                                                     |
| <b>ACCREDITATION (ISO 15189):</b> | 2021-07-08 / 2026-02-02                                                                                                                           |
| <b>TURNAROUND TIME (MAXIMUM):</b> | 3 - 4 months                                                                                                                                      |
| <b>CREATED:</b>                   | 19 Jul 2019 - 12:22                                                                                                                               |
| <b>CHANGED:</b>                   | 01 Mar 2023 - 15:10                                                                                                                               |
| <b>URL:</b>                       | <a href="https://laboboeken.nexuzhealth.com/pboek/internet/GHB%20CME/14756">https://laboboeken.nexuzhealth.com/pboek/internet/GHB%20CME/14756</a> |

-----  
Source URL: [http://gentest.healthdata.be/genetic\\_test/91](http://gentest.healthdata.be/genetic_test/91)

## RELATED CONTENT

---

### Related Diseases

- [Bleeding disorder in hemophilia A carriers](#)
- [Hemophilia A](#)
- [Mild hemophilia A](#)
- [Moderate hemophilia A](#)
- [Severe hemophilia A](#)

### Related Laboratories

- [Centrum Menselijke Erfelijkheid - KUL](#)

### Related Analytes

- [coagulation factor VIII](#)

---

Source URL: [http://gentest.healthdata.be/genetic\\_test/91](http://gentest.healthdata.be/genetic_test/91)