

**GENETIC TEST:**  
**Fanconi anemia (FANCC) (hot spot mutation - c.456+4A>T)**

<b>FULL NAME:</b>	Fanconi anemia (FANCC) (hot spot mutation - c.456+4A>T)
<b>TEST TYPE:</b>	Clinical
<b>TEST SPECIALTY:</b>	Molecular Genetics
<b>TEST PURPOSE:</b>	Carrier diagnosis
<b>SPECIMEN:</b>	Peripheral (whole) blood on EDTA
<b>METHOD CATEGORY:</b>	Targeted variant analysis
<b>METHOD TECHNIQUE:</b>	Single nucleotide extension
<b>RIZIV CODE:</b>	565390-565401
<b>TURNAROUND TIME (MAXIMUM):</b>	1 month
<b>CREATED:</b>	30 Jul 2019 - 09:58
<b>CHANGED:</b>	14 Dec 2022 - 11:59
<b>URL:</b>	<a href="https://labogidsmedgen.uza.be/analyses/fanconi-anemie-type-c">https://labogidsmedgen.uza.be/analyses/fanconi-anemie-type-c</a>

## RELATED CONTENT

---

### Related Diseases

- Fanconi anemia

### Related Laboratories

- Centrum Medische Genetica - UZ Antwerpen

### Related Analytes

- FA complementation group C

---

Source URL: [http://gentest.healthdata.be/genetic\\_test/532](http://gentest.healthdata.be/genetic_test/532)